
DE Gebrauchsanweisung **EN** Instructions for Use **FR** Mode d'emploi
NL Gebruiksaanwijzing **IT** Istruzioni d'uso **TR** Kullanma Kilavuzu
ES Manual de instrucciones

WM 68390 02/2019 DE, EN, FR, NL, IT, TR, ES(MX)








JOYCE SilkGel Full Face

Full Face Mask
vented, non vented

Löwenstein Medical
Technology GmbH + Co. KG
Kronsalsweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.de

€ 0197

1 Bedienung

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1 Operation

- ⇒ Do not close off the exhalation system of the mask.
- ⇒ Only put on the mask for an extended period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who are unable to take it off themselves.
- ⇒ If the mask slips or falls off, the therapy is ineffective.
- Risk of injury if the mask slips!**
- If the mask slips or falls off, the therapy is ineffective.
- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.

2 Introduction

2.1 Intended use

The JOYCE SilkGel Full Face mask is used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device. Use non-vented full face masks only in combinations with therapy devices that have an active exhalation valve.

2.2 Contraindications

The mask may not be used, or may be used only with particular caution, if any of the following symptoms are present:

1. Headgear
2. Forehead cushion
3. Forehead support
4. Retaining ring
5. Elbow
6. Anti-asphyxia valve (only vented)
7. Pressure measurement port seal

12 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

Mask part	Material
Headgear clip	POM (polyoxymethylene)
Forehead support, rotating sleeve	PA (polyamide)
	EL (elasthane), PES (polyester), PU (polyurethane), PA (polyamide), CO (cotton)
Headgear	

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

11 Warranty

Löwenstein Medical gives the customer a limited manufacturer warranty on new original Löwenstein Medical products and any replacement part fitted by Löwenstein Medical in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions are available on the website of the manufacturer. We will also send you the warranty conditions on request.

Please bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine replacement parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

2.3 Side effects

Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

3 Safety

Risk of injury from the supply of oxygen!

Oxygen can become deposited in clothing, bedding and hair. In conjunction with smoking, naked flames or electrical equipment, it can cause fires and explosions.

- ⇒ Do not smoke.
- ⇒ Avoid naked flames.

Risk of injury from re-inhaling CO₂!

If the mask is used incorrectly, CO₂ may be re-inhaled.

4.1 Overview

A diagram of the individual parts can be found on the title page.

4.2 Compatible devices

With many combinations of device, actual pressure in the mask will not correspond to prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that actual pressure in the mask corresponds to therapy pressure.

4.3 Exhalation system

2. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	x	
When washing mask parts (except for mask cushion and forehead cushion), clean thoroughly with a cloth or a soft brush.	x	
Wash headgear by hand.		x

3. Rinse all parts with clear water.

4. Allow all parts to air-dry.
5. Perform a visual inspection.
6. If necessary: replace damaged parts.
7. Re-assemble mask (see Figure 5).

4.4 Anti-asphyxia valve

⚠ WARNING

Risk of asphyxia if anti-asphyxia valve not working properly!

Residues may cause the valve to stick and lead to CO₂ being re-inhaled.

⇒ Check before every use that the openings of the anti-asphyxia valve are clear.

If the therapy device fails, the anti-asphyxia valve opens so that the patient can breathe ambient air.

4.5 Pressure measurement port

The mask has a pressure measurement port, which can be used for measuring pressure or supplying oxygen. If you are not using the connection, close it with the seal for pressure measurement to reach therapy pressure.

5 Hygiene treatment

⚠ WARNING

Risk of injury from inadequate cleaning!

Residues may congest the mask, impair the integrated exhalation system and jeopardize therapy success.

⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

5.1 Clean mask

1. Dismantle mask (see Figure 4).

Fault	Cause	Remedy
Switching pressure Anti-asphyxia valve (only vented)	Mask not correctly adjusted.	Re-adjust mask (see Figure 2).
1 hPa 2 hPa	Mask cushion damaged.	Replace mask cushion.
Quoted two-figure noise emission value to ISO 4871: - Sound pressure level - Sound power level - Uncertainty factor	Patient circuit leaking.	Check connectors, check tubes properly located.
15 dB(A) 23 dB(A) 3 dB(A)	Therapy pressure is not reached. Air escapes at pressure measurement port (only vented).	Close pressure measurement port with plug (included in scope of supply).
up to 12 months ¹	Anti-asphyxia valve incorrectly fitted (only vented).	Fit anti-asphyxia valve correctly (see Figure 5).
EN ISO 17510-2:2009	Anti-asphyxia valve defective (only vented).	Replace anti-asphyxia valve.

- ¹ The mask materials deteriorate if exposed to e.g. aggressive detergents. In individual cases it may be necessary to replace mask parts sooner.
The right to make design modifications is reserved.

9 Pressure/flow curve

The characteristic pressure/flow curve shows the leakage flow as a function of therapy pressure.

8 Technical specifications

Product class to directive 93/42/EEC	IIa
Dimensions (W x H x D)	97 mm x 145 mm x 90 mm
Size S	97 mm x 155 mm x 93 mm
Size M	97 mm x 167 mm x 94 mm
Size L	
Weight	
Size S	159 g
Size M	166 g
Size L	178 g
Therapy pressure range	4 hPa - 40 hPa
Tube connection: tapered connection to EN ISO 5356-1	Ø 22 mm (male) Ø 22 mm (female)
Ventilated non-vented	
Pressure measurement port	Ø 4 mm
Temperature range:	+5 °C to +40 °C
Operation	-20 °C to +70 °C
Storage	

10 Materials

Mask part	Material
Mask cushion, forehead cushion	SI (Silicone), SI (silicone gel)
Plug for pressure measurement port	SI (Silicone)
Retaining ring, anti-asphyxia valve	TPE (thermoplastic elastomer), PP (polypropylene)
Mask body, elbow	PC (polycarbonate), PA (polyamide)

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